Flow cytometry sample preparation

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1. Cell preparation

- 1) 12 well trypsin 100ul 사용하여 세포를 모운다(1300rpm 3 min) (including neutralization step).
- 2) Add 200ul of FACS buffer for washing (V bottom) and resuspension.
- 3) Centrifugation (2000rpm, 3min)
- 4) Discard supernatant.

2. FcγR blocking

- 1) Add 50ul (per well, V bottom) of FcγR blocking solution (1 : 500 = FcγR antibody : FACS buffer) and resuspension.
- 2) Incubation in ice for 30 min
- 3) Add 200ul of FACS buffer (total: 250ul) and centrifugation (2000rpm, 3min).
- 4) Discard supernatant.

3. Flow cytometry staining

*Fixation (with fixation buffer)

- 1) **Fix buffer** 200ul 분주 후 suspension
- 2) Incubation in ice for 20 min
- 3) 20분 이후 2000rpm 3min
- 4) Discard supernatant

*Surface marker staining (dark condition)

- 5) Add 50ul of antibody cocktail (per well, V bottom) (commonly 200-500-fold dilution in FACS buffer) and resuspension.
- 6) Incubation in ice for 30 min.
- 7) Add 200ul of FACS buffer (total: 250ul) and centrifugation (2000rpm, 3min).
- 8) Discard supernatant.

*Permeabilization (with permeabilization buffer)

- 9) Add 200ul of **Perm buffer** and suspension.
- 10) Incubation in ice for 30min
- 11) 30분 이후 2000rpm 3min centrifuge

12) Discard supernatant

***Intracellular marker staining (dark condition)**

13) Add 50ul of antibody cocktail (per well, V bottom) (commonly 200-500-fold dilution in **Perm buffer**) and resuspension.

14) Incubation in ice for 30 min.

15) Add 200ul of **Perm buffer** (total: 250ul) and centrifugation (2000rpm, 3min).

16) Discard supernatant.

17) Add 200ul of FACS buffer and resuspension.

18) FACS tube에 옮기기

4. Flow cytometry operating

1) [Cytometry setting] → Fluorescence selection

2) [Compensation settings] → create compensation controls