

# THP-1 cell culture and differentiation

20210617 Quan HL

## ➤ Materials


- Complete growth medium: RPMI-1640 Medium + 10 % FBS + 1% PS (Penicillin - Streptomycin) + 0.05 mM 2-mercaptoethanol (ATCC에서 권장, 추가하지 않아도 됨)
- T-75 flasks

## ➤ Process


### Step 1: Subculturing

- Resuspension:  $4 \times 10^6$  cells 를 15mL complete growth media 에 띄워 T-75 flasks 에 넣음
- Subculture: 세포 농도가  $\approx 1.2 \times 10^7$  cells 까지 키움 (세포 농도가  $1 \times 10^6$  cells/mL 를 초과하지 않게 함)
- Population doubling time (PDT): 약 26-28 hrs
- Medium Renewal: 3 days
- Culture Conditions: Atmosphere: air 95%; CO<sub>2</sub> 5%
- Temperature: 37.0 °C

DAY 0



DAY 3



### Step 2: Macrophage differentiation

- M0: PMA 100 ng/mL 24hr 분화 후, PMA 30 ng/mL 로 변경하여 48hr 분화 시킴
- M1: PMA 100 ng/mL 24hr 분화 후 PMA 30 ng/mL + LPS 20 ng/mL + hIFN $\gamma$  20 ng/ml 로 변경하여 48hr 분화 시킴
- M2: PMA 100 ng/mL 24hr 분화 후 PMA 30 ng/mL + hIL-4 20 ng/mL 로 변경하여 48hr 분화 시킴